


## P7201-ARMU/P7201-SBMU



①	Seat ass'y	⑤	Controller
	32200134(Pan Seat Ass'y,Black/Gray 24W*20D*25H)		33703023 (DYNAMIC SHARK CONTROLLER)
	32228003 (Pan Seat Ass'y,Black/Gray 24W*20D*25H)		64810005 (JOYSTICK ONLY, DY SHARK,DK-REMA03)
	43701072(PAD, ARMREST(LONG),BLACK)		64210005 (Controller Only,Dy Sharkdk-Pma05).
	43701073(PAD,ARMREST(SHORT),BLACK)	⑥	Shroud ass'y
②	Drive Wheel Ass'y.		30510820 (SHROUD ASS'Y., CANDYAPPLE RED, P710)
	32600423(12" Rear Wheel Ass'y.,4.00-5 Alu. Rim, Black,Pt)		30515928 (SHROUD ASS'Y., BLUE,P710)
③	Caster Wheel with Fork Ass'y.	⑦	Footplate ass'y.
	33100049(Castor Wheel, 2.80/2.50-4 Pufoam Filled Wheel)		33303059 (FOOTPLATE ASS'Y, P720)
④	Motor and gearbox complete, LH and RH.	⑧	Charger unit
	62010030 (Motor/Gearbox Ass'y.A9Y1X04372, P7101, LH)		65102001(Charger, Off Board,Hp8204-6a-Off,115/230VAC)
	62020030 (Motor/Gearbox Ass'y.A9Y1X04372, P7101, RH)		63600008(CHARGER POWERCORD W/PLUG,US SPEC)